

A SMALL SOLUTION FOR INSTANT RESULTS

Navina™ Mini can ease your bowel issues. An intuitive solution to aid incomplete emptying and soiling.



Navina™ Mini

NAVINA MINI CAN SET YOU FREE

Bowel problems can have a significant impact on your health and quality of life, causing distress and embarrassment which can result in avoiding social situations due to the fear of accidents.

These issues are more common than you might think, but unfortunately, they are not commonly discussed. Thankfully, there are products available to help.

One of them is Navina Mini - an easy to use solution that can help with your symptoms:

- ✓ The feeling of not emptying completely
- ✓ Soiling

A SMALL SOLUTION FOR INSTANT RESULTS

Navina Mini is a discreet and intuitive low-volume irrigation solution that can help with symptoms such as incomplete emptying or soiling. With a gentle flush of water, it clears the lower part of the rectum, providing immediate relief and preventing leakage.



A SMART DESIGN KEEPS IT EASY

Shaped to fit your hand
for ergonomic use.

Squeeze-friendly bulb
that holds 185 ml of water.

Rounded tip
for comfortable insertion.

Hygienic
for convenient use.


Extension tube
can be connected for extra reach.

Durable design


- Water container and optional extension tube can be used up to 60 times.
- Rectal cone is single use.
- Navina Mini is 100% free from PVC and phthalates.




HOW TO USE



- 1 Sit down or stand over the toilet. Insert the cone into the rectum.



- 2 Gently flush by squeezing the water container.



- 3 You can use the extension tube if you want to.



For full instructions, scan the code with your smartphone.

Navina Bowel Care. A complete range of bowel care products.

We know bowel care is not a one size fits all, which is why we take pride in being able to offer diversified solutions to meet each individual's needs, all under one brand that promises you the same quality, ease of use and flexibility regardless of product choice.

To see our full range of bowel care products, visit our website by scanning the QR code:



Navina Mini

Intended use

Navina Mini is intended for trans-anal irrigation by instilling a low volume of water into rectum through a single use rectal cone.



Manufacturer:
Wellspect HealthCare,
Aminogatan 1, P.O. Box 14,
SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00.
www.wellspect.com

Indications

Navina Mini is indicated for adults, adolescents and children from 3 years who suffer from incomplete bowel emptying and/or fecal incontinence.



Contraindications

- Known anal or colorectal stenosis
- Active inflammatory bowel disease
- Acute diverticulitis and diverticular abscess
- Colorectal cancer
- Ischemic colitis
- You are within three months of anal or colorectal surgery

As the list may not be exhaustive, healthcare professionals will always consider individual user factors as well. Before use see Instructions for Use.

Reference	Type	Specification
6931000	Navina Mini Start Set	1 water container, 5 rectal cones, 1 extension tube
6930700	Navina Mini Set	1 water container, 30 rectal cones
6930800	Navina Mini Extension Set	1 water container, 30 rectal cones, 1 extension tube
6930900	Navina Mini Refill Set	30 rectal cones

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimize the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcome in a sustainable way, now and for the future.

Wellspect. A Real Difference.

For more information about our products and services, please visit [Wellspect.com](https://www.wellspect.com). Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)

